



Verity®

THE PERFECT SYSTEM TO FIT YOUR RESEARCH

PURIFICATION SYSTEMS

Purification System Designed for your Application

VERITY systems offer a wide array of components to meet all of your purification needs.

VERITY Systems Offer:

- Classic HPLC to flash to CPC, for small to large molecules
- Milligram to kilogram purifications
- Batch to continuous injection techniques
- UV/Vis, ELSD and MS detection options




VERITY® HPLC Purification System with Isocratic Pump and UV/Vis Detector



VERITY® PLC Purification System with Centrifugal Partition Chromatography (CPC) Column and MS Detector

Get the Accuracy and Precision You Need for Trusted Results



TRILUTION® LC Control Software

- **Ease of use** — Systems controlled by intuitive software allowing you to focus on the science
- **Flexibility** — Numerous choices and options let you adapt quickly to your changing research
- **Built-to-last** — Durability backed by Gilson's long history in chromatography

Popular VERITY Configurations

- **VERITY PLC Purification Systems** — Easily switch from preparative to flash chromatography with a simple column change
- **VERITY PLC Systems with CPC Column** — Ideal for natural product purification, PLC with CPC offers improved yield and zero sample loss at less cost
- **VERITY PLC-MS Systems** — Reduce the number of fractions collected and analyzed for additional time and money savings
- **VERITY Oligo Purification System** — Isolate, analyze, and desalt reagent-grade oligonucleotides in one continuous run
- **VERITY GPC Clean-up System** — Improve LC and GC column life and decrease downtime with manual or automated post-extraction clean-up